# Installation Manual

RadiForce®  
RX650  
RX850

Color LCD monitor

## Important

Please read this Installation Manual, and “Instructions for Use” (separate volume) carefully to familiarize yourself with safe and effective usage.

- Please refer to the “Instructions for Use” for basic information ranging from connection of the monitor to a PC to using the monitor.
- The latest Installation Manual is available for download from our web site:

<http://www.eizoglobal.com>

This product has been adjusted specifically for use in the region to which it was originally shipped. If operated outside this region, the product may not perform as stated in the specifications.

---

No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, or otherwise, without the prior written permission of EIZO Corporation. EIZO Corporation is under no obligation to hold any submitted material or information confidential unless prior arrangements are made pursuant to EIZO Corporation's receipt of said information. Although every effort has been made to ensure that this manual provides up-to-date information, please note that EIZO monitor specifications are subject to change without notice.

---

# CONTENTS


<b>CONTENTS</b> .....	<b>3</b>	4-3. Acquiring Life Prediction Data .....	20
<b>Chapter 1 Basic Operations and Functions</b> ....	<b>4</b>	4-4. Setting On/Off for Warning and QC History .....	21
1-1. Basic Adjustment menu operations .....	4	<b>Chapter 5 Power Saving Functions</b> .....	<b>22</b>
1-2. Functions .....	6	5-1. Setting Power Save .....	22
<b>Chapter 2 Adjusting Screens</b> .....	<b>7</b>	5-2. Adjusting the Brightness of Indicator .....	22
2-1. Compatible Resolutions .....	7	5-3. Setting the power saving mode when a person leaves the monitor .....	23
2-2. Setting the Resolution .....	7	<b>Chapter 6 Troubleshooting</b> .....	<b>24</b>
● Monitor setting.....	7	6-1. No picture .....	24
● PC setting.....	7	6-2. Imaging problems .....	25
2-3. Adjusting Color .....	9	6-3. Other problems .....	26
● To select the display mode (CAL Switch mode) .....	9	<b>Chapter 7 Reference</b> .....	<b>28</b>
● To perform advanced adjustments .....	10	7-1. How to Attach the Arm .....	28
● Adjustable settings in each mode .....	10	7-2. How to Install the Wall Mounting Device... ..	29
● To adjust brightness .....	10	7-3. Connecting Multiple PCs .....	30
● To adjust contrast .....	11	● To switch among input signals .....	30
● To adjust color temperature .....	11	● To set input signal selection .....	31
● To adjust gamma .....	12	7-4. Making Use of USB (Universal Serial Bus).....	32
● To adjust hue .....	12	● Required system environment .....	32
● To adjust saturation .....	12	● Connection procedure (Setup of USB function).....	32
● To adjust gain .....	13	7-5. Displaying Monitor Information .....	33
● To set the mode in the monochrome display area .....	13	● Displaying the signal information .....	33
● To set the detection type in the monochrome display area .....	14	● Displaying monitor information.....	33
● To set the highlight area in the monochrome display area .....	14	7-6. Main Default Settings .....	34
<b>Chapter 3 Setting Monitor</b> .....	<b>15</b>	<b>Appendix</b> .....	<b>35</b>
3-1. Skipping Unused Display Modes .....	15	LIMITED WARRANTY (RX650) .....	35
3-2. Setting Language.....	15	LIMITED WARRANTY (RX850) .....	43
3-3. Setting Adjustment Menu Display .....	15	Recycling Information .....	51
3-4. Locking Operation Switches .....	16		
3-5. Enabling/Disabling DDC Communication... ..	16		
3-6. Enabling/Disabling DC5V Output.....	17		
3-7. Showing and hiding the Grayscale Display Error .....	17		
3-8. Return To Default Setting.....	18		
● To reset the color adjustments .....	18		
● To reset all settings .....	18		
<b>Chapter 4 Controlling Monitor Quality</b> .....	<b>19</b>		
4-1. Executing Calibration.....	19		
● Executing SelfCalibration.....	19		
● Viewing SelfCalibration result .....	19		
4-2. Checking Grayscale.....	20		
● Executing Grayscale Check.....	20		
● Viewing Grayscale Check result .....	20		

# Chapter 1 Basic Operations and Functions


## 1-1. Basic Adjustment menu operations

### 1. Displaying Adjustment Menu


1. Touch one of the switches (excluding ).

2. The guide is displayed.


3. Select **| MENU |** to display the Adjustment menu.


### 2. Adjusting/Setting


1. Choose a menu to adjust/set with **| ▲ |** or **| ▼ |**, and then choose **| ENTER |**.

2. Choose an item to adjust/set with **| ▲ |** or **| ▼ |**, and then choose **| ENTER |**.

3. Adjust/set the selected item with **| ▲ |** or **| ▼ |**, and then choose **| OK |**.

### 3. Exiting

1. Choose **| RETURN |** a few times. The Adjustment menu finishes.

# 1-2. Functions

The following table shows all the Adjustment menu's adjustment and setting menus.

Main menu	Item	Reference	
Color*1 	Brightness	"2-3. Adjusting Color" (page 9)	
	Contrast		
	Temperature		
	Gamma		
	Advanced Settings		Hue
			Saturation
			Gain
	Hybrid Area		
	Hybrid Assignment		
Highlight Area	"3-8. Return To Default Setting" (page 18)		
Color Reset			
RadiCS SelfQC 	SelfCalibration	"4-1. Executing Calibration" (page 19)	
	Grayscale Check	"4-2. Checking Grayscale" (page 20)	
	LEA	"4-3. Acquiring Life Prediction Data" (page 20)	
	Settings	"4-4. Setting On/Off for Warning and QC History" (page 21)	
PowerManager 	Power Save	"5-1. Setting Power Save" (page 22)	
	Indicator	"5-2. Adjusting the Brightness of Indicator" (page 22)	
	Presence Sensor	"5-3. Setting the power saving mode when a person leaves the monitor" (page 23)	
Menu Settings 	Language	"3-2. Setting Language" (page 15)	
	Menu Position	"3-3. Setting Adjustment Menu Display" (page 15)	
Tools 	Input	"To switch among input signals" (page 30)	
	Input Selection	"To set input signal selection" (page 31)	
	Mode Skip	"3-1. Skipping Unused Display Modes" (page 15)	
	Signal Info	"7-5. Displaying Monitor Information" (page 33)	
	Monitor Info		
	All Reset	"3-8. Return To Default Setting" (page 18)	
Optional Settings	Key Lock	"3-4. Locking Operation Switches" (page 16)	
	DDC	"3-5. Enabling/Disabling DDC Communication" (page 16)	
	DC5V Output	"3-6. Enabling/Disabling DC5V Output" (page 17)	
	Resolution	"2-2. Setting the Resolution" (page 7)	
	Grayscale Warning	"3-7. Showing and hiding the Grayscale Display Error" (page 17)	

\*1 The adjustment/settings available on "Color" depend on the selected mode. (See "Adjustable settings in each mode" (page 10).)

# Chapter 2 Adjusting Screens

## 2-1. Compatible Resolutions

For details on compatible resolutions, refer to “Compatible Resolutions” in the “Instructions for Use”.

## 2-2. Setting the Resolution

When you connect the monitor to the PC and find that the resolution is improper, or when you want to change the resolution, follow the procedure below.

### ● Monitor setting

#### Procedure

1. Touch  $\odot$  to turn off the monitor.
2. While holding down the  $\text{MODE}$  switch at the far left, touch  $\odot$  for at least two seconds to turn on the monitor.  
The “Optional Settings” menu appears.
3. Choose “Resolution” from “Optional Settings”, and choose  $\text{ENTER}$ .
4. Select one of the following settings using  $\blacktriangle$  or  $\blacktriangledown$ , and then choose  $\text{OK}$ .

	Setting	Details	
RX650	Full	DisplayPort	1640×2048×2
		DVI	1640×2048×2
	6MP	DisplayPort	1536×2048×2
		DVI	1536×2048×2
RX850	8MP	DisplayPort	2048×2160×2
		DVI	2048×2160×2
	QFHD	DisplayPort	1920×2160×2
		DVI	1920×2160×2

5. Select “Finish” with  $\blacktriangle$  or  $\blacktriangledown$ .
6. Choose  $\text{ENTER}$ .  
The “Optional Settings” menu is closed.

### ● PC setting

#### Windows 8 / Windows 7

1. For Windows 8, click the “Desktop” tile on the Start Screen to display the desktop.
2. Right-click the mouse anywhere on the desktop except for icons.
3. From the displayed menu, click “Screen resolution”.
4. On the “Screen Resolution” dialog box, select the monitor.
5. Click “Resolution” to select the desired resolution.
6. Click the “OK” button.
7. When a confirmation dialog box is displayed, click “Keep changes”.

## Windows Vista

1. Right-click the mouse anywhere on the desktop except for icons.
2. From the displayed menu, click "Personalize".
3. On the "Personalization" window, click "Display Settings".
4. On the "Display Settings" dialog box, select the "Monitor" tab and select desired resolution in the "Resolution" field.
5. Click the "OK" button.
6. When a confirmation dialog box is displayed, click "Yes".

## Windows XP

1. Right-click the mouse anywhere on the desktop except for icons.
2. From the displayed menu, click "Properties".
3. When the "Display Properties" dialog box is displayed, click the "Settings" tab and select desired resolution for "Screen resolution" under "Display".
4. Click the "OK" button to close the dialog box.

## Mac OS X

1. Select "System Preferences" from the Apple menu.
2. When the "System Preferences" dialog box is displayed, click "Displays" for "Hardware".
3. On the displayed dialog box, select the "Display" tab and select desired resolution in the "Resolutions" field.
4. Your selection will be reflected immediately. When you are satisfied with the selected resolution, close the window.

## 2-3. Adjusting Color


### ● To select the display mode (CAL Switch mode)

This function allows easy selection of a desired mode according to monitor application.

#### Color Mode

√: Calibration is available

Mode	Purpose	
1-DICOM	Available to display in DICOM mode.	√
2-Custom	Select one of these user modes to define a desired set of parameter settings.	√
3-CAL1	Displays the screen adjusted by calibration software.	√
4-CAL2		
5-CAL3		
6-Hybrid-γ	Automatically identifies the display area for the monochrome and color images on the same screen, and displays each set status.	-
7-sRGB	Suitable for color matching with sRGB compatible peripherals.	-
8-Text	Suitable for displaying texts for word processing or spreadsheets.	-

#### Procedure

1. Touch one of the switches (excluding  $\text{⏻}$  ).  
The guide is displayed.
2. Choose | **MODE** |.  
The mode menu appears at the bottom left of the right screen.
3. The mode among the list is highlighted in turn each time is | **MODE** | chosen.  
You can switch the mode with |  $\blacktriangle$  | or |  $\blacktriangledown$  | while the mode menu is displayed.

→ 1-DICOM ← → 2-Custom ← → 3-CAL1 ← → 4-CAL2 ← → 5-CAL3 ← → 6-Hybrid-γ ← → 7-sRGB ← → 8-Text ← →

#### Note

- The Adjustment menu and the Mode menu cannot be displayed at the same time.

## ● To perform advanced adjustments

“Color” of the Adjustment menu allows you to configure and save a separate set of color settings for each mode.

### Attention

- The monitor will take about 30 minutes before its internal electric components stabilize. Wait 30 minutes or more after turning the monitor power on before starting adjustment.
- The same image may be seen in different colors on multiple monitors due to their monitor-specific characteristics. Make fine color adjustment visually when matching colors on multiple monitors.


### Note

- The values shown in “%” or “K” are available only as reference.

## ● Adjustable settings in each mode

The adjustable settings depend on the mode. (You cannot select any non-adjustable or non-settable function.)

√: Adjustable -: Not adjustable

Icon	Function	CAL Switch mode					
		1-DICOM*1	2-Custom*1	3-CAL1*1 4-CAL2*1 5-CAL3*1	6-Hybrid-y	7-sRGB	8-Text
	Brightness	√	√	√	-	√	√
	Contrast	-	√	-	-	-	√
	Temperature	-	√	-	-	√	√
	Gamma	-	√	-	-	-	√
	Hue	-	√	-	-	-	√
	Saturation	-	√	-	-	-	√
	Gain	-	√	-	-	-	-
	Hybrid Area	-	-	-	√	-	-
	Hybrid Assignment	-	-	-	√	-	-
	Highlight Area	-	-	-	√	-	-
	Color Reset	√	√	√	√	√	√

\*1 If you adjust the calibration in this mode using the calibration kit exclusive for to this monitor (see RadiCS LE /"Accessories" in the "instructions for Use".), only “Brightness” and “Color Reset” are adjustable.

## ● To adjust brightness

The screen brightness is adjusted by changing the brightness of the backlight (light source from the LCD back panel).

### Adjustable range

0 to 100%

### Procedure

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Brightness” from “Color”, and choose | **ENTER** |.
3. Use | **▲** | or | **▼** | for adjustment.
4. Choose | **OK** | to exit.

### Note

- If you feel that the image is too dark even if the brightness is set to 100%, adjust the contrast.

## ● To adjust contrast

The luminance of the screen is adjusted by varying the video signal level.

### **Adjustable range**

0 to 100%

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Contrast” from “Color”, and choose | **ENTER** |.  
“Contrast” appears.
3. Use | ▲ | or | ▼ | for adjustment.
4. Choose | **OK** | to exit.

---

#### **Note**

- In the contrast of 100%, every color gradation is displayed.
  - When adjusting the monitor, it is recommended to perform the brightness adjustment which may not lose the gradation characteristics, prior to the contrast adjustment.
  - Perform the contrast adjustment in the following cases.
    - When you feel the image is too bright even if the brightness is set to 0% (Set the contrast to lower than 100%.)
- 

## ● To adjust color temperature

The color temperature can be adjusted.

The color temperature is normally used to express the hue of “White” and/or “Black” by a numerical value. The value is expressed in degrees “K” (Kelvin).

The screen becomes reddish at a low color temperature, and bluish at a high color temperature, like the flame temperature. The gain preset values are set for each color temperature setting value.

### **Adjustable range**

Native, 6000K-15000K (specified every 100K), User, Standard

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Temperature” from “Color”, and choose | **ENTER** |.
3. Use | ▲ | or | ▼ | for adjustment.
4. Choose | **OK** | to exit.

---

#### **Note**

- “Gain” allows you to perform more advanced adjustment (see [“To adjust gain” \(page 13\)](#)).
  - If you set “Native”, the image is displayed in the preset color of the monitor (Gain: 100% for each RGB).
  - When gain is changed, the color temperature is changed to “User”.
-

## ● To adjust gamma

Gamma can be adjusted. The brightness of the monitor varies depending on the input signal, however, the variation rate is not proportional to the input signal. To keep the balance between the input signal and the brightness of the monitor is called as “Gamma correction”.

### **Adjustable range**

1.6 to 2.7

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Gamma” from “Color”, and choose | **ENTER** |.
3. Use | **▲** | or | **▼** | for adjustment.
4. Choose | **OK** | to exit.

## ● To adjust hue

This function allows you to adjust the hue.

### **Adjustable range**

-100 to 100

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Advanced Settings” from “Color”, and choose | **ENTER** |.
3. Choose “Hue”, and choose | **ENTER** |.
4. Use | **▲** | or | **▼** | for adjustment.
5. Choose | **OK** | to exit.

---

#### **Attention**

- Using this function may make some color gradations unavailable for display.
- 

## ● To adjust saturation

This function allows you to adjust color saturation.

### **Adjustable range**

-100 to 100

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Advanced Settings” from “Color”, and choose | **ENTER** |.
3. Choose “Saturation”, and choose | **ENTER** |.
4. Use | **▲** | or | **▼** | for adjustment.
5. Choose | **OK** | to exit.

---

#### **Attention**

- Using this function may make some color gradations unavailable for display.
- 

#### **Note**

- The minimum value (-100) changes the screen to monochrome.
-

## ● To adjust gain

The brightness of each color component red, green, and blue is called Gain. You can change the hue of “white” by adjusting gain.

### **Adjustable range**

0 to 100%

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Advanced Settings” from “Color”, and choose | **ENTER** |.
3. Choose “Gain”, and choose | **ENTER** |.
4. Choose the color for adjustment from “Red”, “Green”, and “Blue”, and choose | **ENTER** |.
5. Use | **▲** | or | **▼** | for adjustment.
6. Choose | **OK** | to exit.

---

### **Attention**

- Using this function may make some color gradations unavailable for display.

### **Note**

- The value of gain changes with that of color temperature.
  - When gain is changed, the color temperature is changed to “User”.
- 

## ● To set the mode in the monochrome display area

The mode suitable for the monochrome image display area can be set (see [“To select the display mode \(CAL Switch mode\)” \(page 9\)](#)).

### **Adjustable range**

DICOM, CAL1, CAL2, CAL3

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Hybrid Area” from “Color”, and choose | **ENTER** |.
3. Use | **▲** | or | **▼** | for adjustment.
4. Choose | **OK** | to exit.

---

### **Attention**

- You cannot select any mode for which Mode Preset is “Off” (see [“3-1. Skipping Unused Display Modes” \(page 15\)](#)).
- If all “Mode Preset” for DICOM/CAL1/CAL2/CAL3 are set to “Off”, the mode selected last is set.

### **Note**

- The mode excluding the monochrome image display area is “Custom” mode.
  - The Brightness is the same as the value for the mode selected. To change the brightness, use the mode selected for adjustment.
-

## ● To set the detection type in the monochrome display area

You can set the detection type in the monochrome image display area.

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Hybrid Assignment” from “Color”, and choose | **ENTER** |.
3. Select “Full Scan” or “Custom” with | ▲ | or | ▼ |.
4. Choose | **OK** | to exit.

Detection Type	Purpose
Full Scan	Identifies the monochrome image display area automatically.
Custom	Allows you to set the monochrome image display area manually using software.

## ● To set the highlight area in the monochrome display area

You can set On/Off for the highlight area in the monochrome image display area.

### **Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Highlight Area” from “Color”, and choose | **ENTER** |.
3. Select “On” or “Off” with | ▲ | or | ▼ |.
4. Choose | **OK** | to exit.

---

### **Note**

- You can set On/Off for the highlight area with | ■■ |.
-

# Chapter 3 Setting Monitor

## 3-1. Skipping Unused Display Modes

This function selects / deselects modes that should be skipped.

Use this function when all the display modes are not available or when keeping the display mode unchanged.

### Procedure

1. Choose “Tools” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Mode Skip” from “Tools”, and choose | **ENTER** |.
3. Use | **▲** | or | **▼** | to select the mode to change, and choose | **ENTER** |.
4. Select “Skip” or “-” with | **▲** | or | **▼** |.
5. Choose | **OK** | to exit.

### Attention

- You cannot set all the modes disable. Set at least one mode to “-”.

## 3-2. Setting Language

This function allows you to select a language for the Adjustment menu and messages.

### Selectable languages

English, German, French, Spanish, Italian, Swedish, Japanese, Simplified Chinese, Traditional Chinese

### Procedure

1. Choose “Menu Settings” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Language” from “Menu Settings”, and choose | **ENTER** |.
3. Select a language with | **▲** | or | **▼** |.
4. Choose | **OK** | to exit.

## 3-3. Setting Adjustment Menu Display

You can change the display position of the Adjustment menu. The Adjustment menu is displayed on the right screen.

### Procedure

1. Choose “Menu Settings” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Menu Position” from “Menu Settings”, and choose | **ENTER** |.
3. Choose the position with | **▲** | or | **▼** |.
4. Choose | **OK** | to exit.

## 3-4. Locking Operation Switches

This function allows you to lock the adjusted/set status to prevent changing them.

### Procedure

1. Touch  $\mathcal{U}$  to turn off the monitor.
2. While holding down the  $\text{MODE}$  switch at the far left, touch  $\mathcal{U}$  for at least two seconds to turn on the monitor.

The “Optional Settings” menu appears.

3. Choose “Key Lock” from “Optional Settings”, and choose  $\text{ENTER}$ .
4. Choose “Off”, “Menu”, or “All” with  $\blacktriangle$  or  $\blacktriangledown$ , and choose  $\text{ENTER}$ .

Settings	Switches that can be locked
Off (Default setting)	None (All switches are enabled)
Menu	$\text{MENU}$ switches (adjustments/settings using the Adjustment menu)
All	All switches excluding $\mathcal{U}$

5. Select “Finish” with  $\blacktriangle$  or  $\blacktriangledown$ .
6. Choose  $\text{ENTER}$ .

The “Optional Settings” menu is closed.

## 3-5. Enabling/Disabling DDC Communication

This function allows you to enable/disable the DDC communication. Normally, turn it on during use. (This is different from the DDC/CI communication.)

### Attention

- The DDC communication uses signal 1 (DVI) only.

### Procedure

1. Touch  $\mathcal{U}$  to turn off the monitor.
2. While holding down the  $\text{MODE}$  switch at the far left, touch  $\mathcal{U}$  for at least two seconds to turn on the monitor.

The “Optional Settings” menu appears.

3. Choose “DDC” from “Optional Settings”, and choose  $\text{ENTER}$ .
4. Select “On” or “Off” with  $\blacktriangle$  or  $\blacktriangledown$ , and choose  $\text{ENTER}$ .
5. Select “Finish” with  $\blacktriangle$  or  $\blacktriangledown$ .
6. Choose  $\text{ENTER}$ .

The “Optional Settings” menu is closed.

## 3-6. Enabling/Disabling DC5V Output

This function allows you to enable/disable the DC5V output. Normally, turn it off during use.

### Procedure

1. Touch  $\mathcal{U}$  to turn off the monitor.
2. While holding down the | **MODE** | switch at the far left, touch  $\mathcal{U}$  for at least two seconds to turn on the monitor.

The “Optional Settings” menu appears.

3. Choose “DC5V Output” from “Optional Settings”, and choose | **ENTER** |.
4. Select “On” or “Off” with |  $\blacktriangle$  | or |  $\blacktriangledown$  |, and choose | **ENTER** |.
5. Select “Finish” with |  $\blacktriangle$  | or |  $\blacktriangledown$  |.
6. Choose | **ENTER** |.

The “Optional Settings” menu is closed.

## 3-7. Showing and hiding the Grayscale Display Error

If a 6bit (64 color) signal is detected when connecting to the DisplayPort, a grayscale display error can be displayed. If the error message appears, please turn off the main power for the AC adapter, and then turn it on again.

### Attention

- Always use with this set to “ON” (the default setting is “ON”).
- Set this to “OFF” only when using 6bit (64 grayscale) display.

### Procedure

1. Touch  $\mathcal{U}$  to turn off the monitor.
2. While holding down the | **MODE** | switch at the far left, touch  $\mathcal{U}$  for at least two seconds to turn on the monitor.

The “Optional Settings” menu appears.

3. Choose “Grayscale Warning” from “Optional Settings”, and choose | **ENTER** |.
4. Select “On” or “Off” with |  $\blacktriangle$  | or |  $\blacktriangledown$  |, and choose | **ENTER** |.
5. Select “Finish” with |  $\blacktriangle$  | or |  $\blacktriangledown$  |.
6. Choose | **ENTER** |.

The “Optional Settings” menu is closed.

## 3-8. Return To Default Setting

---

There are two different Resets: Color Reset for only resetting the color adjustments to the default settings and All Reset for resetting all settings to the defaults.

---

**Attention**

- After resetting, you cannot undo the operation.

---

**Note**

- Regarding the default setting, see [“7-6. Main Default Settings” \(page 34\)](#).
- 

### ● To reset the color adjustments

This function only resets the color adjustment value for the currently selected mode to the default settings.

**Procedure**

1. Choose “Color” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Color Reset” from “Color”, and choose | **ENTER** |.
3. Select “Execute” with | ▲ | or | ▼ |.
4. Choose | **ENTER** |.

The color adjustments revert to the default settings.

### ● To reset all settings

This function resets all adjustments to the default settings (except the “Input” and “Optional Settings” menus).

**Procedure**

1. Choose “Tools” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “All Reset” from “Tools”, and choose | **ENTER** |.
3. Select “Execute” with | ▲ | or | ▼ |.
4. Choose | **ENTER** |.

All settings excluding the “Input” and “Optional Settings” menus are reset to the default settings.

# Chapter 4 Controlling Monitor Quality

This product has the built-in Integrated Front Sensor, which allows the user to perform periodic calibration and Grayscale Check for the monitor independently.

## Note

- To perform advanced quality control that conforms to medical standards, use the optional monitor quality control software kit “RadiCS UX1”.
- Perform correlation using RadiCS/RadiCS LE to adjust the measurement result of the Integrated Front Sensor to the measurement result of the UX1 sensor for reference. Periodical correlation enables the measurement accuracy of the Integrated Front Sensor equivalent to that of the UX1 sensor. For details, refer to RadiCS/RadiCS LE User’s Manual.
- “RadiCS” allows you to set the adjustment target for the calibration and judgment value for Grayscale Check, and also the schedule for the calibration and Grayscale Check.
- The supplied “RadiCS LE” software also allows you to set the adjustment target and schedule for the calibration. For details, refer to RadiCS LE User’s Manual (on the CD-ROM).

## 4-1. Executing Calibration

This function allows you to execute SelfCalibration and view the latest result. In addition, when the schedule is set using RadiCS/RadiCS LE, the next calibration schedule appears.

### ● Executing SelfCalibration

#### Procedure

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “SelfCalibration” from the “RadiCS SelfQC” menu, and choose | **ENTER** |.
3. Select “Execute” with | **▲** | or | **▼** |.
4. Choose | **OK** | to exit.  
SelfCalibration is executed.

#### Attention

- After SelfCalibration has been executed, perform the Grayscale Check based on the actual temperature and illuminance of the working environment.

### ● Viewing SelfCalibration result

#### Attention

- The result can be viewed when “QC History” is “On”. (See [“4-4. Setting On/Off for Warning and QC History” \(page 21\).](#))

#### Procedure

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “SelfCalibration” from the “RadiCS SelfQC” menu, and choose | **ENTER** |.
3. Select “Result” with | **▲** | or | **▼** |.
4. Select “DICOM”, “Custom”, “CAL1”, “CAL2” or “CAL3” with | **▲** | or | **▼** |.
5. Choose | **OK** | to exit.  
The SelfCalibration result (Max. Error Rate, Last Time, and Error Code) appears.

## 4-2. Checking Grayscale

This function allows you to execute Grayscale Check and view the latest result. In addition, when the schedule is set using RadiCS, the next grayscale check schedule appears.

### ● Executing Grayscale Check

---

**Attention**

- Perform the Grayscale Check based on the actual temperature and illuminance of the working environment.
- 

**Procedure**

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Grayscale Check” from the “RadiCS SelfQC” menu, and choose | **ENTER** |.
3. Select “Execute” with | ▲ | or | ▼ |.
4. Choose | **OK** | to exit.

Grayscale Check is executed.

### ● Viewing Grayscale Check result

---

**Attention**

- The result can be viewed when “QC History” is “On”.
- 

**Procedure**

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Grayscale Check” from the “RadiCS SelfQC” menu, and choose | **ENTER** |.
3. Select “Result” with | ▲ | or | ▼ |.
4. Select “DICOM”, “Custom”, “CAL1”, “CAL2” or “CAL3” with | ▲ | or | ▼ |.
5. Choose | **OK** | to exit.

Up to the five latest Grayscale Check results (Passed/Failed/Canceled/ERR) appear.

## 4-3. Acquiring Life Prediction Data

The data necessary to predict the monitor life can be acquired.

---

**Note**

- Data are collected every 100 hours.
  - When you select “Power Save”, data are collected when the mode switches to Power Save or when the power is shut down after 100 hours.
  - When you select “Routine” or “Power Save”, the Integrated Front Sensor performs measurement during data collection. When you select “Off”, measurement using the Integrated Front Sensor is not performed during data collection.
  - Data are saved after the following times have elapsed.
    - 500, 1000, 2000, 4000, 7000, 10000, 15000, 20000, 25000, 30000 hours
- 

**Procedure**

1. Select “RadiCS SelfQC” from the Adjustment menu, and choose | **ENTER** |.
2. Select “LEA” from “RadiCS SelfQC”, and choose | **ENTER** |.
3. Select “Meas. Frequency” with | ▲ | or | ▼ |.
4. Select “Power Save”, “Routine”, or “Off” with | ▲ | or | ▼ |.
5. Choose | **OK** | to exit.

## 4-4. Setting On/Off for Warning and QC History

This function allows you to set On/Off for Warning and QC History when the Grayscale Check result is Failed.

### **Procedure**

1. Choose “RadiCS SelfQC” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Settings” from the “RadiCS SelfQC” menu, and choose | **ENTER** |.
3. Select “Warning” or “QC History” with | **▲** | or | **▼** |.
4. Select “On” or “Off” with | **▲** | or | **▼** |.
5. Choose | **OK** | to exit.

# Chapter 5 Power Saving Functions

## 5-1. Setting Power Save

This function allows you to set the monitor into the power saving mode according to the PC status. When the monitor has shifted to the power saving mode, the images on the screen are not displayed.

### Attention

- Turning off the AC adapter main power switch or unplugging the power cord completely shuts off power supply to the monitor.
- Devices connected to the USB upstream port and USB downstream port work when the monitor is in power saving mode or when the power switch of the monitor is Off. Therefore, power consumption of the monitor varies with connected devices even in the power saving mode.
- Five seconds before the monitor enters the power saving mode, it displays a preliminary message.

### Procedure

1. Choose “PowerManager” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Power Save” from “PowerManager”, and choose | **ENTER** |.
3. Select “Off”, “High”, or “Low” with | **▲** | or | **▼** |.

Setting	Function
High	Turns off the backlight. The power saving effect is the greatest.
Low	Turns down the backlight to the lowest setting. The time until the operation stabilizes after exiting the Power Save mode next time can be reduced.

4. Choose | **OK** | to exit.

### Power Saving System

#### At DVI signal input

This monitor complies with the “DVI DMPM” standard.

#### At DisplayPort signal input

This monitor complies with the “DisplayPort Standard V1.1a”.

The monitor enters the power saving mode after five seconds in connection with the PC setting.

PC	Monitor	Power Indicator
On	Operating	Green
Power saving mode	Power saving mode	Orange

## 5-2. Adjusting the Brightness of Indicator

The brightness of the power indicator (green) when the screen is displayed can be adjusted (default setting is set to light up when power is turned on, and brightness is set to 4).

### Procedure

1. Choose “PowerManager” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Indicator” from “PowerManager”, and choose | **ENTER** |.
3. Select the indicator brightness “Off” or in the range from 1 to 7 with | **▲** | or | **▼** | as desired.
4. Choose | **OK** | to exit.

## 5-3. Setting the power saving mode when a person leaves the monitor

The sensor on the front of the monitor detects the movement of the operator when Presence Sensor is set to “On”. When a person moves away from the monitor, it shifts to the power saving mode automatically and does not display the images on the screen to reduce the power consumption. When a person comes near the monitor again, it recovers from the power saving mode and displays the images. The sensitivity and time until the power saving mode is activated can be set according to the monitor usage environment and movement of the user.

### Note

- At the time of shifting to power saving mode, the message that indicates shifting is displayed.

### Attention

- This setting cannot be used during “RadiCS SelfQC” execution.

### Procedure

1. Select “PowerManager” from the Adjustment menu, and choose | **ENTER** |.
2. Select “Presence Sensor” from “PowerManager”, and choose | **ENTER** |.
3. Select “On” or “Off” with | **▲** | or | **▼** |.
4. Choose | **OK** | to exit.
5. After selecting “On”, set the details.

Item	Setting Range	Description
Time	5, 30 sec 1, 3, 5, 10, 15, 30, 45, 60 min	Specifies the time from when a person leaves the monitor until the “No Presence detected” message appears. The monitor shifts to the power saving mode about 20 seconds after the message is displayed.
Sensitivity	Level 5 Level 4 Level 3 Level 2 Level 1	When the setting level is high (“Level 5” at maximum), even minor movements of a person are detected, causing the monitor to shift to the power saving mode only infrequently. On the other hand, when the setting level is low (“Level 1” at minimum), only major movements of a person are detected, allowing the monitor to shift to the power saving mode frequently. The Presence Sensor function is activated when the sensor on the front of the monitor detects the movement of heat. Therefore, the sensitivity should be changed depending on the typical user’s clothes or the typical room temperature. If the monitor shifts to or recovers from the power saving mode contrary to your expectations, try changing the setting level.
Reset	-	This function resets all adjustments to the default settings.


### Note

- To use the monitor comfortably, increase the sensitivity for high temperatures, which reduces the detection sensitivity, and reduce the sensitivity for low temperatures, which increases the detection sensitivity.


# Chapter 6 Troubleshooting

If a problem still remains after applying the suggested remedies, contact your local EIZO representative.

## 6-1. No picture

Problem	Possible cause and remedy
<b>1. No picture</b> <ul style="list-style-type: none"> <li>Power indicator does not light.</li> </ul>	<ul style="list-style-type: none"> <li>Check whether the power cord and the adapter power cable are connected properly.</li> <li>Turn the AC adapter main power switch on.</li> <li>Touch .</li> <li>Turn off the AC adapter main power, and then turn it on again a few minutes later.</li> </ul>
<ul style="list-style-type: none"> <li>Power indicator is lighting green.</li> </ul>	<ul style="list-style-type: none"> <li>Increase "Brightness", "Contrast", and/or "Gain" in the Adjustment menu (see <a href="#">"To perform advanced adjustments"</a> (page 10)).</li> <li>Turn off the AC adapter main power, and then turn it on again.</li> </ul>
<ul style="list-style-type: none"> <li>Power indicator is lighting orange.</li> </ul>	<ul style="list-style-type: none"> <li>Switch the input signal (see <a href="#">"To switch among input signals"</a> (page 30)).</li> <li>Move the mouse or press any key on the keyboard.</li> <li>Check whether the PC is turned on.</li> <li>When Presence Sensor is set to "On", the monitor may have shifted to power saving mode. Come near to the monitor.</li> <li>Turn off the AC adapter main power, and then turn it on again.</li> </ul>
<ul style="list-style-type: none"> <li>Power indicator flashes orange and green.</li> </ul>	<ul style="list-style-type: none"> <li>The device that is connected using DisplayPort have a problem. Solve the problem, turn off the monitor, and then turn it on again. Refer to the user's manual of the device that uses DisplayPort connection for details.</li> </ul>
<b>2. The message below appears.</b>	This message appears when the signal is not input correctly even when the monitor functions properly.
<ul style="list-style-type: none"> <li>This message appears when no signal is input. Example:</li> </ul> 	<ul style="list-style-type: none"> <li>The message shown left may appear, because some PCs do not output the signal soon after power-on.</li> <li>Check whether the PC is turned on.</li> <li>Check whether the signal cable is connected properly.</li> <li>Switch the input signal (see <a href="#">"To switch among input signals"</a> (page 30)).</li> </ul>

Problem	Possible cause and remedy
<ul style="list-style-type: none"> <li>The message shows that the input signal is out of the specified frequency range. (Such signal frequency is displayed in magenta.) Example:</li> </ul>  <p>fD: Dot clock fH: Horizontal scan frequency fV: Vertical scan frequency</p>	<ul style="list-style-type: none"> <li>Check whether the PC is configured to meet the resolution and vertical scan frequency requirements of the monitor (see “Compatible Resolutions” in the “Instructions for Use”).</li> <li>Reboot the PC.</li> <li>Select an appropriate setting using the graphics board’s utility. Refer to the manual of the graphics board for details.</li> </ul>

## 6-2. Imaging problems

Problem	Possible cause and remedy
1. The screen is too bright or too dark.	<ul style="list-style-type: none"> <li>Use “Brightness” or “Contrast” in the Adjustment menu to adjust it. (The LCD monitor backlight has a limited life span. When the screen becomes dark or begins to flicker, contact your local EIZO representative.)</li> </ul>
2. Characters are blurred.	<ul style="list-style-type: none"> <li>Check whether the PC is configured to meet the resolution and vertical scan frequency requirements of the monitor (see “Compatible Resolutions” in the “Instructions for Use”).</li> </ul>
3. Afterimages appear.	<ul style="list-style-type: none"> <li>Afterimages are particular to LCD monitors. Avoid displaying the same image for a long time.</li> <li>Use the screen saver or power saving function to avoid displaying the same image for extended periods of time.</li> </ul>
4. Green/red/blue/white dots or defective dots remain on the screen.	<ul style="list-style-type: none"> <li>This is due to LCD panel characteristics and is not a failure.</li> </ul>
5. Interference patterns or pressure marks remain on the screen.	<ul style="list-style-type: none"> <li>Leave the monitor with a white or black screen. The symptom may disappear.</li> </ul>
6. The left and right screens appear misaligned when two screens (left and right) are displayed as one screen.	<ul style="list-style-type: none"> <li>The display timing may be off depending on the graphic card you are using.</li> </ul>

## 6-3. Other problems

Problem	Possible cause and remedy
1. The adjustment menu does not appear.	<ul style="list-style-type: none"> <li>• Check whether the operation lock function works (see “3-4. Locking Operation Switches” (page 16)).</li> </ul>
2. No mode menu is displayed.	<ul style="list-style-type: none"> <li>• Check whether the operation lock function works (see “3-4. Locking Operation Switches” (page 16)).</li> </ul>
3. The operation switches do not work.	<ul style="list-style-type: none"> <li>• Check whether the operation lock function works (see “3-4. Locking Operation Switches” (page 16)).</li> <li>• Check whether there are any water droplets or foreign objects on the switch surface. Gently wipe the switch surface and try operating the switches again with dry hands.</li> <li>• Are you wearing gloves? If so, remove the gloves and try operating the switches again with dry hands.</li> </ul>
<b>4. (When using Presence Sensor)</b> <b>The images remain to be displayed when you move away from the monitor. / The images do not reappear when you come near it.</b>	<ul style="list-style-type: none"> <li>• Check the setting environment of the monitor. Presence Sensor may not work correctly in the following environments. <ul style="list-style-type: none"> <li>- The monitor is located in a place affected by wind.</li> <li>- There is some heat emitting device near the monitor.</li> <li>- There is some obstacle in front of the sensor.</li> </ul> </li> <li>• Check whether any dirt on the sensor. Clean the sensor with a soft cloth.</li> <li>• Check whether you are in front of the monitor. The Presence Sensor function is activated when the sensor on the front of the monitor detects the movement of heat.</li> <li>• In the high temperature conditions, the monitor may not recover from the power saving mode. If the screen is not displayed even after moving the mouse or pressing any key on the keyboard, turn the monitor off and on using its power button. When the screen is displayed, increase the detection sensitivity value (see “5-3. Setting the power saving mode when a person leaves the monitor” (page 23)).</li> </ul>
<b>5. The monitor connected with the USB cable is not detected. / USB devices connected to the monitor do not work.</b>	<ul style="list-style-type: none"> <li>• Check whether the USB cable is connected correctly (see “7-4. Making Use of USB (Universal Serial Bus)” (page 32)).</li> <li>• Change the USB port to another one. If the PC or peripheral devices works correctly by changing the USB port, contact your local EIZO representative. (Refer to the manual of the PC for details.)</li> <li>• Reboot the PC.</li> <li>• If the peripheral devices work correctly when the PC and peripheral devices are connected directly, please contact your local EIZO representative.</li> <li>• Check whether the PC and OS are USB compliant. (For USB compliance of the respective devices, consult their manufacturers.)</li> <li>• Check the PC’s BIOS setting for USB when using Windows. (Refer to the manual of the PC for details.)</li> </ul>

Problem	Possible cause and remedy
<p><b>6. The power turns off during operation and no image is displayed.</b></p>	<ul style="list-style-type: none"> <li>• When the internal temperature is extremely high, the monitor and AC adapter turn off automatically. If the power switch is turned on again while the temperature remains high, the power turns off again in one minute.</li> <li>• Change the installation location or environment, and then turn it on again a few minutes later. <ul style="list-style-type: none"> <li>- The ventilation slot must not be closed.</li> <li>- There is no high temperature object around.</li> </ul> </li> <li>• If the power is not turned on even after the installation location or environment has changed, contact your local EIZO representative.</li> </ul> <p><b>Attention</b></p> <ul style="list-style-type: none"> <li>• Do not block the ventilation slots for the monitor or AC adapter, and do not use the monitor or AC adapter in a poorly-ventilated area.</li> </ul>
<p><b>7. The operation sound of the fan for the AC adapter is noticeable.</b></p>	<ul style="list-style-type: none"> <li>• The AC adapter is equipped with a fan to avoid rise in the internal temperature. The operation sound may be heard depending on the installation location. However, it is not a problem.</li> <li>• Change the installation location of the AC adapter or the orientation of the fan ventilation slot.</li> </ul>
<p><b>8. The Hybrid Gamma function cannot be used.</b></p>	<ul style="list-style-type: none"> <li>• Check whether the input resolutions of the left and right screens are the same.</li> <li>• This function can be used only when the input resolutions of the left and right screens are the same or the signals are input to either screen only.</li> </ul>
<p><b>9. SelfCalibration/Grayscale Check failure</b></p>	<ul style="list-style-type: none"> <li>• Refer to the error code table in the “Instructions for Use”.</li> <li>• If an error code that does not appear in the error code table is displayed, contact your local EIZO representative.</li> </ul> <p><b>Attention</b></p> <ul style="list-style-type: none"> <li>• Do not touch the Integrated Front Sensor.</li> </ul>
<p><b>10. Time is not displayed correctly in the history of SelfCalibration/Grayscale Check.</b></p>	<ul style="list-style-type: none"> <li>• Execute Monitor Detection using RadiCS/RadiCS LE. For details, refer to RadiCS/RadiCS LE User’s Manual.</li> </ul>

# Chapter 7 Reference

## 7-1. How to Attach the Arm

Another manufacturer's arm (or another manufacturer's stand) can be attached by removing the stand section of the monitor.

### Attention

- When attaching an arm or stand, follow the instructions of their User's Manual.
- When using another manufacturer's arm or stand, confirm the following in advance and select one conforming to the VESA standard.
  - Clearance between the screw holes: 100mm × 100mm
  - Strong enough to support weight of the monitor unit (except the stand) and attachments such as cables.
- Use the screws as described below.
  - Screws fixing the stand to the monitor
- When using an arm or stand, attach it to meet the following tilt angles of the monitor.
  - Up 45°, down 45°
- Connect the cables after attaching an arm or stand.
- Do not move the removed stand up and down. Doing so may result in injury or device damage.
- The monitor and arm or stand are heavy. Dropping them may result in injury or equipment damage.

### Procedure


**1. Lay the LCD monitor on a soft cloth spread over on a stable surface with the panel surface facing down.**

**2. Remove the stand.**

Prepare a screwdriver. Loosen the screws (four) securing the unit and the stand with the screwdriver.

**3. Attach the arm or stand to the monitor.**

Attach it using the screws removed in Step 2.

## 7-2. How to Install the Wall Mounting Device

The monitor and AC adapter can be attached to a wall mounting device.

### Attention

- Ask an installation contractor for the attachment.
- Prepare a wall mounting device separately.
- When installing a wall mounting device, follow the instructions of its User's Manual.
- Use the screws as described below.
  - Monitor: Screws fixing the stand to the monitor
  - AC adapter: Screws on the backside (3 pcs.)
- Attach the AC adapter in the correct orientation in accordance with the procedure.
- After the attachment of the AC adapter, be sure to tighten the screws on the backside.
- Connect the cables after attaching an arm or stand.
- Do not move the removed stand up and down. Doing so may result in injury or device damage.
- The monitor, AC adapter, and wall mounting device are heavy. Dropping them may result in injury or equipment damage.

### Procedure

**1. Lay the LCD monitor on a soft cloth spread over on a stable surface with the panel surface facing down.**

**2. Remove the stand.**


Prepare a screwdriver. Loosen the screws (four) securing the unit and the stand with the screwdriver.

**3. Install the wall mounting device and attach the monitor to it.**

Attach it using the screws removed in Step 2.

**4. Attach the AC adapter to the wall mounting device.**

Remove the three screws from the backside of the AC adapter.


Attach the AC adapter as shown below using the screws removed. Do not attach the AC adapter in the wrong orientation.

## 7-3. Connecting Multiple PCs

The product has multiple connections to PCs and allows you to switch to one of the connections for display.

### Connection examples


### ● To switch among input signals

#### Procedure

1. Choose "Tools" from the Adjustment menu, and choose **| ENTER |**.
2. Choose "Input" from "Tools", and choose **| ENTER |**.
3. Select "DVI - DVI", "DVI - DP", "DP - DP", or "DP - DVI" with **| ▲ |** or **| ▼ |**.
4. Choose **| ENTER |** to exit.

When the signal is switched, the active input port name appears at the top right corner of the screen.

## ● To set input signal selection

Setting	Function
Auto	The monitor automatically identifies the connectors carrying computer input signals and displays the respective images on the left and right screens. ( Example: DisplayPort 1 - DisplayPort 2 / DisplayPort 1 - DVI 2 / DVI 1 - DisplayPort 2 / DVI 1 - DVI 2 ) When a PC is turned off or enters the power saving mode, the monitor automatically displays another signal.
Manual	The monitor detects only the PC's signals currently displaying automatically. Select an active input signal with "Input".

### **Procedure**

1. Choose "Tools" from the Adjustment menu, and choose | **ENTER** |.
2. Choose "Input Selection" from "Tools", and choose | **ENTER** |.
3. Select "Auto" or "Manual" with | ▲ | or | ▼ |.
4. Choose | **ENTER** | to exit.

---

### **Note**

- When "Auto" is selected for "Input Selection", the monitor's power saving function works only when all PCs are in the power saving mode.
-

## 7-4. Making Use of USB (Universal Serial Bus)

This monitor has a hub compatible with USB. Connected to a PC compatible with USB or another USB hub, this monitor functions as a USB hub allowing connection to peripheral USB devices.

### ● Required system environment


- A PC equipped with a USB port or another USB hub connected to a USB compatible PC
- Windows 2000/XP/Vista/7/8 or Mac OS 9.2.2 and Mac OS X 10.2 or later
- USB cable

#### Attention

- This monitor may not work depending on PC, OS or peripheral devices to be used. For USB compatibility of peripheral devices, contact their manufactures.
- When the monitor is in power saving mode, devices connected to the USB downstream port will still work. Therefore, power consumption of the monitor varies with connected devices even in the power saving mode.
- When the AC adapter main power switch is Off, device connected to the USB port will not operate.

### ● Connection procedure (Setup of USB function)

1. Connect the monitor first to a PC using the signal cable, and run the PC.
2. Connect the supplied USB cable between the USB downstream port of a USB compatible PC (or another USB hub) and the monitor's USB upstream port.  
The USB hub function is set up automatically upon connection of the USB cable.
3. Connect the peripheral USB device to the USB downstream port of the monitor.

## 7-5. Displaying Monitor Information

### ● Displaying the signal information

This function displays the information about the current input signals displayed.

#### Procedure

1. Choose “Tools” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Signal Info” from “Tools”, and choose | **ENTER** |.

The “Signal Info” appears.

(Example)

```
Tools > Signal Info
Signal-1 DVI      Signal-2 DP
1280 X 1024      1280 X 1024
fD: 108.0MHz     fD: 108.0MHz
fH: 63.9kHz      fH: 63.9kHz
fV: 60.0Hz       fV: 60.0Hz
```

### ● Displaying monitor information

This function displays the information about the monitor.

#### Procedure

1. Choose “Tools” from the Adjustment menu, and choose | **ENTER** |.
2. Choose “Monitor Info” from “Tools”, and choose | **ENTER** |.

The “Monitor Info” appears.

This function displays the model name, serial number, firmware version, and usage time.

(Example)

```
Tools > Monitor Info
RadiForce RX650
S/N: 00000001
Version: 10000-10000-00000
Usage Time
Product:      1h
Backlight:    1h
```

---

#### Attention

- The usage time is not always “0” when you purchase the monitor due to factory inspection.
-

## 7-6. Main Default Settings

**CAL Switch Mode : The default display mode setting is 1-DICOM.**

### RX650

Mode	Brightness	Temperature	Gamma value
1-DICOM	400 cd/m <sup>2</sup>	7,500K	DICOM
2-Custom	Approx. 400 cd/m <sup>2</sup>	7,500K	2.2
3-CAL1	300 cd/m <sup>2</sup>	7,500K	DICOM
4-CAL2	300 cd/m <sup>2</sup>	7,500K	2.2
5-CAL3	400 cd/m <sup>2</sup>	6,500K	DICOM
6-sRGB	Approx. 250 cd/m <sup>2</sup>	6,500K	Standard
7-Text	Approx. 100 cd/m <sup>2</sup>	6,500K	2.2

### RX850

Mode	Brightness	Temperature	Gamma value
1-DICOM	500 cd/m <sup>2</sup>	7,500K	DICOM
2-Custom	Approx. 500 cd/m <sup>2</sup>	7,500K	2.2
3-CAL1	400 cd/m <sup>2</sup>	7,500K	DICOM
4-CAL2	400 cd/m <sup>2</sup>	7,500K	2.2
5-CAL3	500 cd/m <sup>2</sup>	6,500K	DICOM
6-sRGB	Approx. 250 cd/m <sup>2</sup>	6,500K	Standard
7-Text	Approx. 100 cd/m <sup>2</sup>	6,500K	2.2

### Others

Input* <sup>1</sup>	DVI-DVI
Presence Sensor	Off
Power Save	High
Input Selection	Auto
DDC* <sup>1</sup>	On
DC5V Output* <sup>1</sup>	Off
Resolution* <sup>1</sup>	RX650: Full RX850: 8MP
Grayscale Warning	On
Menu Position	Center
Language	English

\*<sup>1</sup> These functions cannot be initialized with the "All Reset" function (see [page 18](#)).

# Appendix

## LIMITED WARRANTY (RX650)

### LIMITED WARRANTY

EIZO Corporation (hereinafter referred to as “**EIZO**”) and distributors authorized by EIZO (hereinafter referred to as the “**Distributors**”) warrant, subject to and in accordance with the terms of this limited warranty (hereinafter referred to as the “**Warranty**”), to the original purchaser (hereinafter referred to as the “**Original Purchaser**”) who purchased the product specified in this document (hereinafter referred to as the “**Product**”) from EIZO or Distributors, that EIZO and Distributors shall, at their sole discretion, either repair or replace the Product at no charge if the Original Purchaser becomes aware within the Warranty Period (defined below) that (i) the Product malfunctions or is damaged in the course of normal use of the Product in accordance with the description in the instruction manual attached to the Product (hereinafter referred to as the “**User’s Manual**”), or (ii) the LCD panel and brightness of the Product cannot maintain the recommended brightness specified in the User’s Manual in the course of normal use of the Product in accordance with the description of the User’s Manual.

The period of this Warranty is five (5) years from the date of purchase of the Product (hereinafter referred to as the “**Warranty Period**”), provided that the brightness of the Product shall be warranted only if the Product has been used within the recommended brightness described in the User’s Manual and the Warranty Period for the brightness is limited to five (5) years from the date of purchase of the Product subject to the usage time being less than or equal to 20,000 hours (the brightness is 400 cd/m<sup>2</sup> and color temperature is 7500K) or less than or equal to 30,000 hours (the brightness is 300 cd/m<sup>2</sup> and color temperature is 7500K). EIZO and Distributors shall bear no liability or obligation with regard to the Product in relation to the Original Purchaser or any third parties other than as provided under this Warranty.

EIZO and Distributors will cease to hold or store any parts of the Product upon expiration of seven (7) years after the production of such parts is discontinued. In repairing the monitor, EIZO and Distributors will use renewal parts which comply with our QC standards.

The Warranty is valid only in the countries or territories where the Distributors are located. The Warranty does not restrict any legal rights of the Original Purchaser.

Notwithstanding any other provision of this Warranty, EIZO and Distributors shall have no obligation under this Warranty whatsoever in any of the cases as set forth below:

- (a) Any defect of the Product caused by freight damage, modification, alteration, abuse, misuse, accident, incorrect installation, disaster, faulty maintenance and/or improper repair by third party other than EIZO and Distributors;
- (b) Any incompatibility of the Product due to possible technical innovations and/or regulations;
- (c) Any deterioration of the sensor;
- (d) Any deterioration of display performance caused by the deterioration of expendable parts such as the LCD panel and/or backlight, etc. (e.g. changes in color, changes in color uniformity, defects in pixels including burnt pixels, etc.);
- (e) Any deterioration of the Product caused by the use at higher brightness than the recommended brightness described in the User’s Manual;
- (f) Any defect of the Product caused by external equipment;
- (g) Any defect of the Product on which the original serial number has been altered or removed;
- (h) Any normal deterioration of the product, particularly that of consumables, accessories, and/or attachments (e.g. buttons, rotating parts, cables, User’s Manual, etc.); and
- (i) Any deformation, discoloration, and/or warp of the exterior of the product including that of the surface of the LCD panel.

To obtain service under the Warranty, the Original Purchaser must deliver the Product, freight prepaid, in its original package or other adequate package affording an equal degree of protection, assuming the risk of damage and/or loss in transit, to the local Distributor. The Original Purchaser must present proof of purchase of the Product and the date of such purchase when requesting services under the Warranty.

The Warranty Period for any replaced and/or repaired product under this Warranty shall expire at the end of the original Warranty Period.

EIZO OR DISTRIBUTORS ARE NOT RESPONSIBLE FOR ANY DAMAGE TO, OR LOSS OF, DATA OR OTHER INFORMATION STORED IN ANY MEDIA OR ANY PART OF ANY PRODUCT RETURNED TO EIZO OR DISTRIBUTORS FOR REPAIR.

EIZO AND DISTRIBUTORS MAKE NO FURTHER WARRANTIES, EXPRESSED OR IMPLIED, WITH RESPECT TO THE PRODUCT AND ITS QUALITY, PERFORMANCE, MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE. IN NO EVENT SHALL EIZO OR DISTRIBUTORS BE LIABLE FOR ANY INCIDENTAL, INDIRECT, SPECIAL, CONSEQUENTIAL OR OTHER DAMAGE WHATSOEVER (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFIT, BUSINESS INTERRUPTION, LOSS OF BUSINESS INFORMATION, OR ANY OTHER PECUNIARY LOSS) ARISING OUT OF THE USE OR INABILITY TO USE THE PRODUCT OR IN ANY CONNECTION WITH THE PRODUCT, WHETHER BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE, EVEN IF EIZO OR DISTRIBUTORS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. THIS EXCLUSION ALSO INCLUDES ANY LIABILITY WHICH MAY ARISE OUT OF THIRD PARTY CLAIMS AGAINST THE ORIGINAL PURCHASER. THE ESSENCE OF THIS PROVISION IS TO LIMIT THE POTENTIAL LIABILITY OF EIZO AND DISTRIBUTORS ARISING OUT OF THIS LIMITED WARRANTY AND/OR SALES.

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# LIMITED WARRANTY (RX850)

## LIMITED WARRANTY

EIZO Corporation (hereinafter referred to as “EIZO”) and distributors authorized by EIZO (hereinafter referred to as the “Distributors”) warrant, subject to and in accordance with the terms of this limited warranty (hereinafter referred to as the “Warranty”), to the original purchaser (hereinafter referred to as the “Original Purchaser”) who purchased the product specified in this document (hereinafter referred to as the “Product”) from EIZO or Distributors, that EIZO and Distributors shall, at their sole discretion, either repair or replace the Product at no charge if the Original Purchaser becomes aware within the Warranty Period (defined below) that (i) the Product malfunctions or is damaged in the course of normal use of the Product in accordance with the description in the instruction manual attached to the Product (hereinafter referred to as the “User’s Manual”), or (ii) the LCD panel and brightness of the Product cannot maintain the recommended brightness specified in the User’s Manual in the course of normal use of the Product in accordance with the description of the User’s Manual.

The period of this Warranty is five (5) years from the date of purchase of the Product (hereinafter referred to as the “Warranty Period”), provided that the brightness of the Product shall be warranted only if the Product has been used within the recommended brightness described in the User’s Manual and the Warranty Period for the brightness is limited to five (5) years from the date of purchase of the Product subject to the usage time being less than or equal to 20,000 hours (the brightness is 500cd/m<sup>2</sup> and color temperature is 7500K) or less than or equal to 30,000 hours (the brightness is 400cd/m<sup>2</sup> and color temperature is 7500K). EIZO and Distributors shall bear no liability or obligation with regard to the Product in relation to the Original Purchaser or any third parties other than as provided under this Warranty.

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- (a) Any defect of the Product caused by freight damage, modification, alteration, abuse, misuse, accident, incorrect installation, disaster, faulty maintenance and/or improper repair by third party other than EIZO and Distributors;
- (b) Any incompatibility of the Product due to possible technical innovations and/or regulations;
- (c) Any deterioration of the sensor;
- (d) Any deterioration of display performance caused by the deterioration of expendable parts such as the LCD panel and/or backlight, etc. (e.g. changes in brightness uniformity, changes in color, changes in color uniformity, defects in pixels including burnt pixels, etc.);
- (e) Any deterioration of the Product caused by the use at higher brightness than the recommended brightness described in the User’s Manual;
- (f) Any defect of the Product caused by external equipment;
- (g) Any defect of the Product on which the original serial number has been altered or removed;
- (h) Any normal deterioration of the product, particularly that of consumables, accessories, and/or attachments (e.g. buttons, rotating parts, cables, User’s Manual, etc.); and
- (i) Any deformation, discoloration, and/or warp of the exterior of the product including that of the surface of the LCD panel.

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---

## Recycling Information

This product, when disposed of, is supposed to be collected and recycled according to your country's legislation to reduce environmental burden. When you dispose of this product, please contact a distributor or an affiliate in your country.

The contact addressees are listed on the EIZO website below.

<http://www.eizoglobal.com>

For recycling information for customers in Switzerland, please refer to the following website.

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